# CAT

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**MCS 77 / MCS 78**

# Table of Contents

<b>1.</b>	<b>User Instructions</b>	<b>5</b>
1.1	Important Instructions for your safety	5
1.2	Danger symbols in this operating manual	6
1.2.1	Danger symbols	6
1.2.2	Danger levels	6
1.2.3	Warning signs on the device	7
<b>2.</b>	<b>General safety warnings and instructions</b>	<b>7</b>
<b>3.</b>	<b>Scope of Delivery</b>	<b>10</b>
<b>4.</b>	<b>General Information</b>	<b>10</b>
<b>5.</b>	<b>Intended Use</b>	<b>10</b>
<b>6.</b>	<b>Setting up and Starting up the Instrument</b>	<b>11</b>
6.1	Unpacking the Instrument	11
6.2	Setting up the Instrument	11
6.3	Electrical Connection	11
6.4	Rear Connections	12
6.4.1	Mains Connection	12
6.4.2	Connection of the Pt100 Probe	13
6.4.3	Connection of the RS485 Interface	14
6.5	Probe Stand	14
<b>7.</b>	<b>Operation of the device</b>	<b>15</b>
7.1	Description of the Controls	15
7.2	Switching the Device On and Off	16
<b>8.</b>	<b>Heating</b>	<b>17</b>
8.1	Temperature Control with External Pt100 Probe	17
8.1.1	Setting up the Heating Volume	18
8.1.2	Hotplate temperature Limit	19
8.2	Temperature Control without External Pt100 Probe (Hotplate Temperature Control)	20
8.3	Ramp Function	20
<b>9.</b>	<b>Stirring</b>	<b>21</b>
<b>10.</b>	<b>Timed Operation</b>	<b>21</b>
10.1	Timer Function	21
10.2	Multitimer Function	22
<b>11.</b>	<b>Safety Functions</b>	<b>23</b>
11.1	Safety Stir Function	23
11.2	Safety Temperature	24
11.2.1	Safety Temperature with External Pt100 Probe	24
11.2.2	Safety Temperature without External Pt100 Probe	25
11.3	Differential Alarm	25
11.4	Out of Liquid Check	26
11.5	Failure of External Pt100 Probe	26
11.6	Hotplate Temperature Limit	27
11.7	Self-Monitoring of the Instrument	27
<b>12.</b>	<b>Device Settings</b>	<b>28</b>
12.1	Setup Menu	28
12.2	Settings Menu	29
12.2.1	Slave Address	29
12.2.2	Baud Rate	29
12.2.3	Temperature Unit	29
12.2.4	Differential Alarm	30
12.2.5	Out of Liquid Check	30
12.2.6	Plate Limit	30
12.2.7	Plate Auto Set	30
12.2.8	Safety Auto Set	31
12.2.9	Ask Volume	31
12.2.10	Motor Auto On	31
12.2.11	Safety Stir Time	31
12.3	PID Menu	32
12.3.1	Thermal Resistance	32
12.3.2	Container Type	32
12.4	Calibration Menu	32

12.4.1	Calibration for Pt100A.....	33
12.4.2	Reset Calibration .....	33
12.5	Default Settings Menu.....	34
12.5.1	Default Settings.....	34
<b>13.</b>	<b>Display Messages .....</b>	<b>34</b>
13.1	Messages at Starting .....	34
13.2	Messages at Shutdown.....	35
13.3	Information and Warnings.....	35
13.4	Off Conditions .....	36
<b>14.</b>	<b>RS485 Interface .....</b>	<b>37</b>
14.1	Settings of the RS485 Interface.....	38
14.2	Pin assignment of the RS485 Connector.....	38
14.3	Connection between PC and Instrument.....	39
14.4	Addressing the Instrument.....	39
14.5	RS485 Data Transfer.....	40
14.5.1	Format of an RS485 Command.....	40
14.5.2	Format of an RS485 Handshake .....	40
14.6	RS485 Commands.....	41
14.7	RS485 Examples .....	45
14.7.1	Example 1 – Setting up a Value .....	45
14.7.2	Example 2 – Reading a Value .....	46
<b>15.</b>	<b>Maintenance and Cleaning .....</b>	<b>47</b>
<b>16.</b>	<b>Warranty and Liability .....</b>	<b>47</b>
<b>17.</b>	<b>Dismantling, Transport and Storage.....</b>	<b>48</b>
17.1	Dismantling .....	48
17.2	Transport and Storage.....	48
<b>18.</b>	<b>Technical Data.....</b>	<b>48</b>
<b>19.</b>	<b>Repair Return Form .....</b>	<b>50</b>

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# 1. User Instructions

## 1.1 Important Instructions for your safety

- Every user must read and understand this manual completely before use. Failure to do so can result in serious injury or death.
- Comply with all safety and accident-prevention regulations applicable to laboratory work.
- Follow general instructions for hazard prevention and general safety instructions, e.g. wear protection clothing, eye protection and gloves.
- This operating manual is part of the product. Thus, it must always be easily accessible.
- This instruction sheet does not purport to address all of the safety problems which might result from the use of this device, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult their authorized safety advisors and establish appropriate health and safety practices and then determine the application of regulatory limitations prior to use.
- Enclose this operating manual when transferring the device to another place.
- If this manual is lost, please request another one. Please contact your dealer or

**Ingenieurbüro CAT  
M. Zipperer GmbH  
Etzenbach 16  
D-79219 Staufen  
Tel. ++49-7636-7803-0  
info@cat-ing.de**

## 1.2 Danger symbols in this operating manual

The safety instructions in this manual appear with the following danger symbols and danger levels:

### 1.2.1 Danger symbols

	<b>Hazard point</b>		<b>Electrical shock</b>
	<b>Risk of fire</b>		<b>Explosion</b>
	<b>Crushing</b>		<b>Hot surface</b>
	<b>Bio hazard</b>		<b>Magnetic field</b>
	<b>Chemical hazard</b>		<b>Material damage</b>

### 1.2.2 Danger levels

<b>▲ DANGER</b>	Will lead to severe injuries or death
<b>▲ WARNING</b>	May lead to severe injuries or death
<b>▲ CAUTION</b>	May lead to light to moderate injuries
<b>NOTICE</b>	May lead to material damage

### 1.2.3 Warning signs on the device

	<b>▲WARNING</b>	This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.
	<b>▲WARNING</b>	Beware of the effect of the <b>magnetic field</b> on cardiac pacemakers and data storage media.
	<b>▲WARNING</b>	This symbol indicates that there are <b>hot surfaces</b> on the device. It refers to the hotplate of the device which can reach up to 500°C during operation causing a hazard resulting in serious burns. <ul style="list-style-type: none"> <li>● Avoid contact with the hotplate, vessel and heating medium.</li> <li>● Avoid storage of flammable substances near the instrument.</li> </ul>

## 2. General safety warnings and instructions

	<b>▲DANGER</b> Risk of explosion.	<ul style="list-style-type: none"> <li>● Do not operate the device in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof.</li> <li>● Do not use this device for processing any substances which could generate an explosive atmosphere.</li> <li>● Do not use this device to process any explosive or highly reactive substances.</li> <li>● Do not use flammable or explosive substances near the instrument.</li> </ul>
	<b>▲DANGER</b> Electric shock as a result of penetration of liquid.	<ul style="list-style-type: none"> <li>● Do not allow any liquids to penetrate the inside of the housing</li> <li>● Switch off the device and disconnect the power plug before starting cleaning or disinfection work. The On/Off Switch on the device does not disconnect the device from the power source.</li> <li>● Use only closed tubes.</li> <li>● Only plug the device back in if it is completely dry, both inside and outside.</li> </ul>
	<b>▲WARNING</b> Risk from incorrect supply voltage	<ul style="list-style-type: none"> <li>● Only connect the device to a AC power source with an protective earth (PE).</li> <li>● Only connect the device to voltage sources which correspondent to the electrical requirements on the type label.</li> </ul>

	<p><b>▲WARNING</b> Electric shock due to damage to device or mains cable</p> <ul style="list-style-type: none"> <li>● Only connect the device to the mains supply if the device and the mains cable are undamaged</li> <li>● Only use devices that have been properly installed or repaired.</li> <li>● In case of danger, disconnect the device from the mains supply by pulling the power plug from the mains socket or by using the isolating device intended for this purpose (e.g. emergency stop switch)</li> </ul>
	<p><b>▲WARNING</b> Lethal voltage inside the device</p> <ul style="list-style-type: none"> <li>● Do not open the device.</li> <li>● Ensure that the housing is always closed and undamaged so that no parts inside the housing can be contacted by accident.</li> <li>● The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.</li> <li>● Do not allow any liquids to penetrate the inside of the housing.</li> <li>● Repairs are only to be carried out by trained service technicians.</li> </ul>
	<p><b>▲WARNING</b> Damage to health due to corrosive or aggressive chemicals</p> <ul style="list-style-type: none"> <li>● Observe all markings on the reagent bottles.</li> <li>● Always check the instrument for leaks and air bubbles. Special attention should be directed to determine that all push-ons, threaded connections and suction tubes are firmly in place before beginning operation. Leaking solutions may endanger persons and materials</li> <li>● Use proper connecting vessels, protective clothing and gloves.</li> <li>● Avoid splashes</li> <li>● When dispensing, maintain a physical distance between the instrument and the body.</li> <li>● Dangerous and fuming chemicals must be dispensed in a fume hood.</li> <li>● Only employ the instrument for the purpose intended by the manufacturer, and particularly within the resistance limits of the instrument. If in doubt, contact your supplier, or the manufacturer's factory representative at the phone number shown at the front page of this operating instruction.</li> <li>● Always use the instrument in such a manner that neither the operator, nor any other person is endangered.</li> </ul>
	<p><b>▲WARNING</b> Damages to health due to infectious liquids and pathogenic germs.</p> <ul style="list-style-type: none"> <li>● When handling infectious liquids and pathogenic germs, observe the national regulations, the biological security level of your laboratory, the material safety data sheets and the manufacturer's application notes.</li> <li>● Wear personal protective equipment</li> <li>● For comprehensive regulations about handling germs or biological material of the risk group II or higher, please refer to the "Laboratory Biosafety Manual" in its respectively current valid version from the World Health Organisation</li> </ul>

	<p><b>▲WARNING</b> <b>Damage to health due to contaminated device and accessories</b></p> <p>In the following cases, sample material can be released:</p> <ul style="list-style-type: none"> <li>- improperly sealed tubes</li> <li>- unstable tubes</li> <li>- high vapour pressure of the content so that the seal of the tubes can spring open</li> <li>- damaged sealings</li> <li>- smashed glass tubes</li> </ul> <ul style="list-style-type: none"> <li>● Only mix in closed tubes</li> <li>● Observe the nationally prescribed safety environment when working with hazardous, toxic and pathogenic samples. Pay particular attention to personal protective equipment (gloves, clothing, goggles, etc.), extraction, and the safety class of the lab.</li> <li>● Decontaminate the device and the accessories before storage and shipping.</li> </ul>
	<p><b>▲WARNING</b> <b>Risk of fire</b></p> <ul style="list-style-type: none"> <li>● Do not use this device to process any highly flammable liquids</li> </ul>
	<p><b>▲CAUTION</b> <b>Poor safety due to inadequate fixing of the unit</b></p> <ul style="list-style-type: none"> <li>● Ensure that the unit is firmly attached to a solid stand.</li> </ul>
	<p><b>▲CAUTION</b> <b>Poor safety due to incorrect accessories and spare parts.</b></p> <p>The use of accessories and spare parts other than recommended by Ingenieurbüro CAT, M. Zipperer GmbH may impair the safety, function and precision of the device. Ingenieurbüro CAT, M. Zipperer GmbH cannot be held liable or accept any liability for damage resulting from the use of incorrect or non-recommended accessories and spare parts, or from the improper use of such equipment.</p> <ul style="list-style-type: none"> <li>● Only use accessories and spare parts recommended by Ingenieurbüro CAT, M. Zipperer GmbH</li> </ul>
	<p><b>▲CAUTION</b> <b>Crush hazard due to moving parts</b></p> <ul style="list-style-type: none"> <li>● Do not replace any consumables as long as the device is running.</li> <li>● Do not open the coverage as long as the device is running</li> </ul>

### 3. Scope of Delivery

Unpack the delivered contents carefully and check to see that nothing appears physically damaged or is missing. Your delivery should contain the following:

- 1 Magnetic Hotplate Stirrer MCS 77 (230V) PN: 60276-50
- or
- 1 Magnetic Hotplate Stirrer MCS 77 (115V) PN: 60276-51
- or
- 1 Magnetic Hotplate Stirrer MCS 78 (230V) PN: 60277-50
- or
- 1 Magnetic Hotplate Stirrer MCS 78 (115V) PN: 60277-51
- 1 Instruction Manual

#### Optional:

- Pt100 Temperature Probe PN: 60278-00
- Probe stand Ø8, 350 mm long, with thread M6, 10 mm long PN: 60663-00
- Probe Clamp PN: 8B005624
- Clamp holder PN: 60668-00
- RS485-RS232 adapter for RS485 communication PN: 61703-00
- RS485 cable for connection between 2 devices
- consisting of: 2x PN: 31112-42
- required length in meters x PN: 10728-66
- RS232 to USB Adapter PN: 30244-01
- 9 pol. Sub-D extension cord PN: 30275-51

### 4. General Information

Our Magnetic Hotplate Stirrers are designed according to protection class 1. They were manufactured and tested according to DIN EN 61010. According to these regulations, the devices are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

### 5. Intended Use

These Magnetic Hotplate Stirrers are mainly to be used to heat and stir hydrous solutions and oil not containing inflammable parts in glass vessels. They are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

## 6. Setting up and Starting up the Instrument

### 6.1 Unpacking the Instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

If the instrument is not damaged and all parts are complete you may start to operate the device after reading the instruction manual.

### 6.2 Setting up the Instrument

Please place the device on a fire-proof horizontal and even surface. A minimum distance of 50cm to inflammable materials should always be observed.


	<b>⚠ DANGER</b> The device may not be operated in explosion-prone areas.
	<b>⚠ WARNING</b> The power cable may not come in contact with the hotplate.
	<b>⚠ WARNING</b> The device is not to be used without supervision.
	<b>⚠ WARNING</b> When connecting the instrument to an AC power outlet ensure that your local supply voltage matches the indication on the instrument.
	<b>⚠ WARNING</b> Beware of the effect of the magnetic field on cardiac pacemakers and data storage media.

### 6.3 Electrical Connection

	<p><b>⚠ WARNING</b> -A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is standard on all instruments. For North America instruments feature standard US plugs (NEMA Pub.No.WDI1961 ASA C 73.1. 1961 page 8 15A 125V).</p> <p>-When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.</p> <p>-The instrument is earthed as supplied. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!</p> <p>-When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.</p>
---	---

## 6.4 Rear Connections

Following connections are located on the rear panel of the device:



### 6.4.1 Mains Connection

Connect the mains power cable to the mains power connector. Ensure that your local supply voltage matches the indication on the instrument.

## 6.4.2 Connection of the Pt100 Probe

The 5-pin DIN Connector on the rear of the instrument is intended for the connection of a Pt100 temperature probe (PN: 60278-00). Pt1000 probes are not suitable.

### Pin assignment of the connector for the Pt100 probe:

	<b>CAUTION</b> Ensure that the cable of the Pt100 probe does not come in contact with the hotplate.
	<b>NOTICE</b> Use Pt100 probes with electrically isolated tips only. Probes that allow for an electrical connection to the metal chassis of the device may distort temperature measurement.
	<b>NOTICE</b> Ensure that the Pt100 probe tip is submersed to a depth of at least 50mm in the heating medium to allow accurate measurement.

### 6.4.3 Connection of the RS485 Interface


On the rear of the device is a connector for a RS485 Interface. For an easy use with a PC, you can use an optional RS485 to RS232 adapter (PN: 61703-50). With this adapter, the device can communicate with a PC.

For more information on this interface please refer to chapter 14-RS485 Interface.

	<b>CAUTION</b> Ensure that the cable does not come in contact with the hotplate.
	<b>NOTICE</b> Only connect original accessories to the RS485 jack! Any other use may result in damage of the device or damage to third party devices!

### 6.5 Probe Stand


Use the M6 thread on the back of the device to mount an optional probe stand (PN: 60663-00). Then it is possible with an additional probe clamp (PN: 8B005624) to clamp the Pt100 probe on a secure and easy way.

## 7. Operation of the device

### 7.1 Description of the Controls

On the control panel, the various device functions can be displayed and set. The different parameters, operating instructions and device states are displayed on the graphic display and the LEDs. With the buttons on the left side you can switch on or off the device as well as the hotplate and motor function. Use the buttons below the display to switch between the different views and set the displayed values.

Bedienelement	Funktion
Display	The display is used for showing the various temperature and speed values as well as different messages and setup menus.
● Heater	This LED indicates that the hotplate is currently being exposed to an energy pulse.
● Hot !	While the hotplate is hotter than 60°C this LED lights up.
● Timer	The Timer LED lights when the timer function is activated. If the multi-timer function is activated, the LED flashes.
● Ramp	The Ramp LED lights when a temperature ramp is set.
	If the device is in standby mode, the device can be switched on by pressing the I/O key. In switched on mode you can press the I/O key again to switch the device back into standby mode. The I/O LED lights up when the device is switched on and flashes in standby mode.
	Press the Plate On/Off key to switch the hotplate on or off. The Plate On/Off LED lights up when the hotplate is turned on.
	Press the Motor On/Off key to switch the motor on or off. The Motor On/Off LED lights up when the hotplate is turned on.
	Press the Plate Temp key to access the plate menu. In this view you can see the actual as well as the set hotplate temperature.
	Press the Probe Temp key to access the probe menu. In this view you can see the actual as well as the set probe temperature. Pressing and holding the key will open the Volume menu.

	Press the Ramp key to access the ramp menu. In this view you can see the actual ramp value.
	Press the Motor Speed key to access the motor menu. In this view you can see the actual as well as the set motor speed.
	Press the Safety Temp key to access the safety menu. In this view you can see the actual safety temperature.
	Press the Timer key to access the timer menu. In this view you can see the actual timer value. Pressing and holding the key will open the Multitimer menu.
	With the Up/+ key the different setpoints can be incremented or navigate in menus upwards.
	With the Down/- key the different setpoints can be decremented or navigate in menus downwards.

## 7.2 Switching the Device On and Off

	<b>▲WARNING</b> The On/Off Switch on the device does not disconnect the device from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.
--	---

Make sure that you have read the instruction manual carefully before operating the device. It is also necessary to check whether the device is in accordance with the setup instructions above.

	If the device is in standby mode, the device can be switched on by pressing the I/O key. In switched on mode you can press the I/O key again to switch the device back in standby mode.
---	---

During switching on and off the device, various messages can appear in the display. These messages give information to the user about different device states. For more information refer to chapter 13-Display Messages.

### Note:

When you switch off the device in standby mode, all values are stored (e.g. the set hotplate temperature). Next time you switch on the device these values are restored.

## 8. Heating

The Magnetic Hotplate Stirrer allows 2 different modes of operation:

- **Temperature control with external Pt100 probe (recommended mode)**
- **Temperature control of the hotplate (no Pt100 probe connected)**

For exact temperature control of a liquid, the use of an external Pt100 probe is recommended to regulate the liquid temperature as accurately as possible.

	<b>▲WARNING</b> Always switch off the device after use!
	<b>▲CAUTION</b> As long as the Heater LED is flashing or lights, the hotplate is supplied with energy. The Hot LED lights when the hotplate is hotter than 60°C. You may get burned!
	<b>▲CAUTION</b> <b>After switching off the device the hotplate stays hot!</b> Therefore, after switching off the device, a warning message appears on the display if the hotplate exceeds 60°C and the power cord is still plugged in.

### 8.1 Temperature Control with External Pt100 Probe

	<b>▲WARNING</b> Please ensure that the cable of the Pt100 probe does not come in contact with the hotplate.
--	---

We recommend the use of an external Pt100 probe (Pt100 Temperature Probe, PN: 60278-00) to allow for accurate temperature control of the heating medium. To heat up a medium, following steps are required:

- Connect the Pt100 probe with the DIN connector on the rear of the Magnetic Hotplate Stirrer
- Place the Pt100 probe in the vessel with the heating medium (min. 5 cm immersed)

	Press the I/O key to switch on the device.
	The device asks about the amount of the medium to be heated. Set the desired value with the Up/+ and the Down/- keys.
	Confirm the set value by pressing the probe key. Another press of the probe key enters the probe menu.
	Set the desired temperature value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.
	After making sure that all parameters are set correctly (also Safety Temp, etc.), the heating process can be started. To do this, press the Plate On/Off key to switch on the hotplate. To switch off the hotplate press the Plate On/Off key again.

**Note:**

When a Pt100 probe is connected, the “Plate Temp” value does not refer to the set temperature of the hotplate, but to the plate limit temperature of the hotplate (see chapter 8.1.2-Hotplate temperature Limit).

**Note:**

When a Pt100 probe is connected to the device, the set temperature of the “Probe Temp” will be set automatically to 0°C and the hotplate will be switched off if it is turned on.

**8.1.1 Setting up the Heating Volume**

For the temperature control of the device, it makes a big difference whether it heats-up a small or a big amount of liquid. So it is necessary to give the device the information about the amount of the heating medium, to achieve the best compromise between temperature accuracy and the time duration of the heating process. The device asks you automatically about the amount of the heating medium (Volume) when it is switched on or a Pt100 probe is connected to the device during operation.

You also can change the amount of the heating medium (Volume) during operation mode manually:

	Press and hold the Probe Temp key to access the medium menu. In this view you can see the actual volume of the liquid.
	Set the desired volume with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.

<p><b>⚠WARNING</b></p>	<p>This value reflects the heat capacity of the liquid on the hotplate. An increase of this value will result in a faster heat up time. If a too high value is entered (much higher than the real existing amount of liquid) this can result in an overshoot and oscillation of the liquid temperature. A too low setting of this value will result in a very slow reaction of the system and a longer than necessary heat up time.</p>
------------------------	---

### 8.1.2 Hotplate temperature Limit

If a Pt100 sensor is connected, the set value "Plate Temp" doesn't refer to the set value of the hotplate, but to the temperature limit of the hotplate. This means that during the heating process the hotplate doesn't rise above the set Plate Limit value.

It should be noted that if the hotplate temperature limit has been limited to a too low value, the liquid may not reach the set temperature (Probe Set) or an extremely long heating time is required.

To avoid this and if no security risk exists, the hotplate temperature limit (Plate Limit) should be set to its max. possible value. So the Magnetic Stirrer has the opportunity to choose the most effective hotplate temperature for a fast heating-up.

Change this value only if there is a security risk (i.e. liquid with a low flash point).

	Press the Plate Temp key to access the plate limit menu. In this view you can see the actual hotplate temperature limit.
	Set the desired plate temperature limit with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.

**Note:**

With every change of the Probe Temperature the hotplate limit is automatically set to its maximum! If necessary, the value must be set back lower.

Should the hotplate temperature be permanently limited, it must be limited as described in chapter 11.7. Hotplate Temperature Limit.

## 8.2 Temperature Control without External Pt100 Probe (Hotplate Temperature Control)

If you want set the temperature of the hotplate directly and not the temperature of the medium, this is only possible when no Pt100 sensor is connected.

	<p>Press the I/O key to switch on the device.</p>
	<p>Press the Plate Temp key to access the Plate menu. In this view you can see the actual as well as the set hotplate temperature.</p>
	<p>Set the desired probe temperature limit with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.</p>
	<p>After making sure that all parameters are set correctly (also Safety Temp, etc.), the heating process can be started. To do this, press the Plate On/Off key to switch on the hotplate. To switch off the hotplate press the Plate On/Off key again.</p>

	<p><b>CAUTION</b> If a Pt100 sensor is connected, the set value Plate Temp doesn't refer to the set value of the hotplate, but to the temperature limit of the hotplate (see chapter 8.1.2-Hotplate temperature Limit).</p>
--	---

### Note:

If the Pt100 sensor is being unplugged, then the set temperature of the hotplate is automatically set to 0°C.

## 8.3 Ramp Function

The device offers the possibility to heat up a liquid (with external Pt100) or the hotplate (without external Pt100) with a defined heating rate (Temperature Ramp). By default this feature is turned off. So the device is trying to reach the set value of the temperature as fast as possible. For example, if the Ramp is set to a value of 10°C/h, the device is trying to reach the set value of the temperature by heating up the liquid or the hotplate with a defined heating rate of 10°C/h.

To turn on the ramp function, the following steps are necessary:

	<p>Press the Ramp key to access the ramp menu. In this view you can see the actual ramp value.</p>
	<p>Set the desired ramp value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.</p>

## 9. Stirring

With this magnetic stirrer, fluids can be stirred in a vessel with a magnetic rod. Make sure that the used vessel is permeable to magnetic field lines (e.g. glass, ceramics or stainless steel).

To adjust the stirring speed and switch on the motor following steps are required:

	Press the Motor Speed key to access the motor menu. In this view you can see the actual as well as the set motor speed value.
	Set the desired motor speed with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.
	Press the Motor On/Off key to switch on the motor. To switch off the motor press the Motor On/Off key again.

	<b>CAUTION</b> Use a suitable vessel (e.g. Erlenmeyer flask) and beware of high stirring speeds and high liquid levels to avoid splashing. Make sure that suitable protective clothing and eyewear is used.
---	---

### Note:

Optimal magnetic coupling is acquired with stirring bars with lengths of 30 – 70 mm and diameters of 5 – 10 mm. The appropriate stirring bar for an application must be determined by the user experimentally.

## 10. Timed Operation

### 10.1 Timer Function

This device has an integrated timer function. The timer function switches off the device after the adjusted time. If the timer expires, the hotplate turns off. The device continues stirring until the adjusted Safety Stir Time expires (see chapter 11.1-Safety Stir Function) and then turns off into standby mode by itself.

The timer function can be switched on by the following way:

	Press the Timer key to access the timer menu. In this view you can see the actual as well as the set timer value.
	Set the desired timer value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the probe temp key.

If you want to cancel the current timer, you have to set the timer value to “Off” as described above.

## 10.2 Multitimer Function

In addition to the standard Timer-Function which turns off the device after a specified time, the device has a Multitimer Function. With this function you can create individual temperature and stirring profiles. The Multitimer Function has 5 programmable steps (Step 1 Set ... Set Step 5). The following parameters can be set for each timer step:

- **Time:** -Execution time (hours : minutes : seconds).  
or  
-Speed reached, Pt100 Probe temperature reached or Plate temperature reached.  
or  
-Step is not executed (Off).
- **Plate:** -Hotplate setpoint in °C or F (without Pt100 Probe).  
or  
-Hotplate temperature limit in °C or F (with Pt100 Probe).
- **Probe:** -Pt100 Probe setpoint in °C or F.
- **Ramp:** -Hotplate temperature Ramp in °C/h or F/h (without Pt100 Probe).  
or  
-Pt100 Probe temperature ramp in °C/h or F/h (with Pt100 Probe).
- **Motor:** -Motor Speed in rpm.

**Note:**

If you select one of the options Speed Reached, Probe Reached or Plate Reached in time setting, the selected step runs until the setpoint is reached. Then it will proceed with the next step.

The Multitimer can also be programmed in the Options menu that a programmed profile runs repeated multiple times or endless:

- **Cycles:** -Endless (The set steps are repeated automatically until the user cancels the operation).  
or  
-Number of repetitions of the steps.

In addition, the behavior of the device after the timer has expired can be set in the Options menu. The following options can be selected:

- **Hold Settings:** -The device continues operation with the settings of the last executed step.
- **Plate Off:** -The device switches off the hotplate. The motor continues operation with the settings of the last executed step.
- **All Off:** -The device switches off the hotplate and motor.
- **Device Off:** -The device will be switched off in standby mode.

To activate the Multitimer Function, the following steps must be performed:

	<p>Press and hold the Timer Temp key to access the Multitimer menu.</p>
	<p>Use the Up/+ and Down/- keys to navigate through the displayed Multitimer menu. Select the menu item --- BACK --- in each submenu to navigate back to main view of the Multitimer menu.</p>
	<p>To change a selected value, press the Timer key. The selected value starts to flash.</p>
	<p>Set the desired value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the timer key.</p>
	<p>After all values have been set and checked, select START TIMER in the main view of the multi-timer menu to activate the function. The motor and the hotplate are automatically turned on. To deactivate the Multitimer Function timer press and hold the Timer key for 3 seconds.</p>

	<p><b>Note:</b> The Multitimer Function can't be combined with the normal Timer function.</p>
--	---

## 11. Safety Functions

The following section describes the different security mechanisms of the device.

### 11.1 Safety Stir Function

The device has different functions (Timer function, ...) and security mechanisms (Safety Temp, ...) that can switch off the device automatically. After such an automatic switch-off, it is often useful to continue stirring the liquid for a certain time to prevent heat accumulation and helping cooling down the liquid. Therefore the device has an integrated Safety Stir Function. This function allows the user to set a desired time period (Safety Stir Time). After an automatic switch-off, the device shuts down the hotplate and continues stirring for the set Safety Stir Time. The Safety Stir Function is only active when the motor is switched on while the heating process. If the Safety Stir Function is active, the display shows an appropriate message with the remaining Safety Stir Time. Once the time has elapsed the device switches off. The user can also switch off the device manually by pressing the I/O key.

To adjust the Safety Stir Time or to switch off this function, please refer to chapter 12.2.11-Safety Stir Time.

	<p><b>NOTICE</b> For serious shutdowns as, for example, a hardware defect of the device, the device is switched off immediately.</p>
---	--

<p><b>Note:</b> If the status of the device allows, the Safety Stir Function continues stirring after every automatic shutdown</p>
--

## 11.2 Safety Temperature

In the case of unexpected heating of the liquid or malfunction of the device, the user can set a Safety Temperature (Safety Temp). If the Safety Temperature is exceeded, the device switches off automatically. If a Pt100 probe is connected to the device, the Safety Temp refers to the temperature of the heating medium, otherwise to the hotplate temperature.

	<b>WARNING</b>	If there is no security risk, the user should set the Safety Temperature to a value as high as possible to avoid an unintentional switch off of the device in normal operation. In safety critical applications the proper use of the safety temperature is advised.
	<b>NOTICE</b>	If a Pt100 probe is connected to the device, the Safety Temp refers to the temperature of the heating medium, otherwise to the hotplate temperature.

### 11.2.1 Safety Temperature with External Pt100 Probe

The device checks if the measured temperature values of the Pt100 probe are lower than the set Safety Temp. When the temperature of the Pt100 probe exceeds the adjusted Safety Temp, the device is switching off the hotplate and continues stirring the adjusted Safety Stir Time as described in chapter 11.1-Safety Stir Function. Afterwards it switches off automatically. If however the sensor temperature still increases to more than 15% above the Safety Temperature, the device will immediately switch off and does not wait until the Safety Stir Time has expired.

If the Probe Temperature is changed by the user, the Safety Temperature is adjusted automatically 15°C above the new set Probe Temperature. You can switch off this automatic function as described in chapter 12.2.8-Safety Auto Set. When the device changes the Safety Temperature automatically, the new value is displayed a short time to inform the user.

To adjust the Safety Temp manually, please follow the steps below:

	Press the Safety Temp key to access the safety menu. In this view you can see the actual Safety Temp value.
	Set the desired value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the Safety temp key.

	<b>WARNING</b>	If the user changes the set temperature of the external Pt100 probe (Probe Temp) above the set safety temperature (Safety Temp), the device automatically adjusts the safety temperature 15°C above the new Probe Temp. When the device changes the Safety Temperature automatically, the new value is displayed.
<p><b>Note:</b> To set the safety temperature (Safety Temp) as close to the actual set value (Probe Temp) as possible, set the probe temperature first, before changing the safety temperature to the desired value manually.</p>		

### 11.2.2 Safety Temperature without External Pt100 Probe

The device checks if the measured temperature values of the hotplate are lower than the set Safety Temp. When the temperature of the hotplate exceeds the adjusted Safety Temp, the device is switching off the hotplate and continues stirring the adjusted Safety Stir Time as described in chapter 11.1-Safety Stir Function. Afterwards it switches off the device automatically. If however the sensor temperature still increases to more than 15% above the Safety Temperature the device will immediately switch off and does not wait until the Safety Stir Time has expired.

If the Plate Temperature is changed by the user, the Safety Temperature is adjusted automatically 15°C above the new set Plate Temperature. You can switch off this automatic function as described in chapter 12.2.8-Safety Auto Set. When the device changes the Safety Temperature automatically, the new value is displayed a short time to inform the user.

To adjust the Safety Temp manually, please follow the steps below:

	Press the Safety Temp key to access the safety menu. In this view you can see the actual Safety Temp value.
	Set the desired value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the Safety temp key.

	<p><b>▲WARNING</b> If the user changes the set temperature of the hotplate (Plate Temp) above the set safety temperature (Safety Temp), the device automatically adjusts the safety temperature 15°C above the new Plate Temp. When the device changes the Safety Temperature automatically, the new value is displayed.</p>
<p><b>Note:</b> To set the safety temperature (Safety Temp) as close to the actual set value (Plate Temp) as possible, set the plate temperature first, before changing the safety temperature to the desired value manually.</p>	

### 11.3 Differential Alarm

This safety function detects an extreme temperature drop of the temperature measured in the medium to be heated (external Pt100 probe). Such a fast temperature drop may occur, for example, if the glass breaks or the Pt100 probe falls out of the liquid. In this case the device is switching off the hotplate and continues stirring the adjusted Safety Stir Time as described in chapter 11.1-Safety Stir Function. Afterwards it switches the device off automatically.

To adjust the sensitivity of the Differential Alarm or switching off this function, please refer to chapter 12.2.4-Differential Alarm.

<p><b>Note:</b> If the device switches off incorrectly, although the Pt100 probe is immersed in the liquid correctly, in certain circumstances it may be necessary to reduce the sensitivity of the Differential Alarm. For example, this could be necessary when the user wants to add a cold liquid into the heated medium.</p>
---

## 11.4 Out of Liquid Check

The Out of Liquid function monitors, whether the Pt100 probe is immersed in the liquid in the vessel on the hotplate. The device checks, whether the Pt100 probe temperature changes in relation to the temperature changes of the hotplate. If the Out of Liquid function detects no increase in the liquid temperature over a certain time, although the hotplate temperature is rising, the device is switching off the hotplate and continues stirring the adjusted Safety Stir Time as described in chapter 11.1-Safety Stir Function. Afterwards it switches the device off automatically.

To adjust the sensitivity of the Out of Liquid function or switching this function off, please refer to chapter 12.2.5-Out of Liquid Check.

**Note:**

If the device is switching off incorrectly, although the Pt100 probe is immersed in the liquid correctly, in certain circumstances it may be necessary to reduce the sensitivity of the Out of Liquid function. For example, this could be necessary when operating with larger volumes of liquid or with liquids that have a very high heat capacity.

**Note:**

This Safety feature is only available, if the set Pt100 probe temperature is a minimum of 10°C higher than the actual liquid temperature.

## 11.5 Failure of External Pt100 Probe

If a Pt100 probe is connected to the device, it monitors if the measured temperatures are in a valid value range. If this is not true or the Pt100 probe is disconnected in operating mode, the device recognizes this. In this case the device is switching off the hotplate and continues stirring the adjusted Safety Stir Time as described in chapter 11.1-Safety Stir Function. Afterwards it switches off automatically.

## 11.6 Hotplate Temperature Limit

The maximum adjustable hotplate temperature can be limited in two different ways:

### Permanent limitation of the hotplate temperature in Setup Menu:

If you want to limit the maximum hotplate temperature permanently, e.g. when the hotplate in a laboratory or a classroom should never rise above 300°C, you have to set this in the Setup Menu as described in section 12.2.6-Plate Limit. If the value in this submenu is limited to 300°C, the user can not set hotplate temperatures over 300°C. The device behaves as if the maximum hotplate temperature is only 300°C. This setting is stored even after switching off and on the device.

### One-time limitation of the hotplate temperature for a single heating process:

The other possibility is limiting the hotplate temperature only for a single heating process. This can be done as described in section 8.1.2-Hotplate temperature Limit. This setting will be reset to its maximum value every time you change the temperature for the heating medium (Probe Temp) and must be re-adjusted if necessary.

#### **▲WARNING**

The use of the hotplate temperature limitation results in a longer heat-up time of the heating medium. Therefore in normal operation, the hotplate temperature limit should be set to the maximum value to achieve optimal heat-up times. Setting the hotplate temperature limit to the maximum value allows the microprocessor to freely select the optimum hotplate temperature to reach the set probe temperature as fast as possible. However, when operating with flammable liquids it may be required to limit the maximum allowed hotplate temperature to a safe value.

## 11.7 Self-Monitoring of the Instrument

The device has additional, internal security functions to monitor and switch off itself automatically, in case of malfunction. For example the device monitors the internal temperature and the internal data communication. When the internal microcontroller hangs up, it will be restarted by the internal watchdog timer. Additionally there are further security mechanisms integrated, to monitor some malfunctions of the heating stage, the temperature probe of the hotplate, as well as wrong calibration values. If a malfunction occurs, the device is switching off the hotplate and continues stirring the adjusted Safety Stir Time as described in chapter 11.1-Safety Stir Function. Afterwards it switches off automatically. For serious shutdowns as, for example, a hardware defect in the device, the device is switched off immediately.

## 12. Device Settings

Besides the two basic functions of heating and stirring, the device offers several additional settings that can be changed via the setting menu.

### 12.1 Setup Menu

The device offers the opportunity to make several basic settings. In the Setup Menu you can select the different submenus, which contains the various setting parameters. The menu can be left via the menu item „BACK“.

- Settings Menu (Section: 12.2-Settings Menu)
- PID Menu (Section: 12.3-PID Menu)
- Calibrate Menu (Section: 12.4-Calibration Menu)
- Default Settings (Section: 12.5-Default Settings Menu)
- BACK

In each submenu different parameters can be selected and changed. The various parameters are described in the corresponding chapters.

To enter the Setup Menu the following steps are required:

	<p>Press the I/O key to switch on the device.</p>
	<p>Directly after pressing and releasing the I/O key, press and hold the Timer key until "Setup Menu" is displayed</p>
	<p>Navigate between the different submenus with the Up/+ and the Down/- keys.</p>
	<p>Press the Timer key to select the shown submenu. The submenu can be left via the menu item „BACK“.</p>

Navigation through the submenus and changing the parameters can be carried out as described below:

	<p>Navigate between the different parameters with the Up/+ and the Down/- keys.</p>
	<p>Press the Timer key to change the selected parameter. The parameter starts to flash.</p>
	<p>Set the desired value with the Up/+ and the Down/- keys. The set value is entered automatically after a short time or can be confirmed manually by pressing the Safety temp key.</p>

## 12.2 Settings Menu

In the Settings Menu you can change various basic settings of the device. This submenu can be reached over the Setup Menu as described in Section 12.1-Setup Menu and contains the following setting parameters:

- Slave Address (Section: 12.2.1-Slave Address)
- Baud Rate (Section: 12.2.2-Baud Rate)
- Temperature Unit (Section: 12.2.3-Temperature Unit)
- Differential Alarm (Section: 12.2.4-Differential Alarm)
- Out of Liquid Check (Section: 12.2.5-Out of Liquid Check)
- Plate Limit (Section: 12.2.6-Plate Limit)
- Plate Temp Auto (Section: 12.2.7-Plate Auto Set)
- Safety Auto Set (Section: 12.2.8-Safety Auto Set)
- Ask Volume (Section: 12.2.9-Ask Volume)
- Motor Auto On (Section: 12.2.10-Motor Auto On)
- Safety Stir Time (Section: 12.2.11-Safety Stir Time)
- BACK

### 12.2.1 Slave Address

If you want to operate several devices on one RS485 bus, each device needs its unique slave address, so that the devices can differentiate between their commands. For more information, see chapter 14.4-Addressing the Instrument.

To adjust the RS485 Slave Address, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Slave Address".

**Adjustment range:** 1 ... 255

**Default value:** 1

### 12.2.2 Baud Rate

With this setting the baud rate of the RS485 port can be changed. It is recommended not to change this value.

To adjust the RS485 baud rate, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Baud Rate".

**Adjustment range:** 1200, 2400, 4800, 9600baud

**Default value:** 9600baud

### 12.2.3 Temperature Unit

The temperature unit can be switched between Celsius and Fahrenheit. It affects both the temperature values on the display, as well as the values read via the RS485 interface.

To adjust the Temperature Unit, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Temp. Unit".

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**Default value:** Celsius

#### 12.2.4 Differential Alarm

This safety function registers an extreme drop in temperature of the measured medium temperature (external Pt100 probe). For more information, see section 11.3-Differential Alarm. To adjust the sensitivity of the Differential Alarm, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Diff. Alarm".

**Adjustment range:** 1% ... 100% (low: 1 ... 39%, normal: 40 ... 69%, high: 70 ... 100%)

**Default value:** 90%

#### 12.2.5 Out of Liquid Check

The Out of Liquid function monitors, whether the Pt100 probe is immersed into the liquid on the hotplate. For more information, see section 11.4-Out of Liquid Check.

To adjust the sensitivity of the Out of Liquid function, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Out Of Liq."

**Adjustment range:** 0 ... 100% (off: 0%, low: 1% ... 39%, normal: 40 ... 69%, high: 70 ... 100%)

**Default value:** 40%

#### 12.2.6 Plate Limit

By changing this value, the maximum adjustable hotplate temperature can be limited. For more information, see section 11.7. Hotplate Temperature Limit

To adjust the maximum allowed hotplate temperature, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Plate Limit".

**Adjustment range:** MCS 77: 50 ... 330°C (Off: 330°C)

MCS 78: 50 ... 440°C (Off: 440°C)

**Default value:** MCS 77: 330°C (Off)

MCS 78: 440°C (Off)

#### 12.2.7 Plate Auto Set

If a Pt100 probe is connected to the device and the temperature of the heating medium (Probe Temp) is changed, the maximum hotplate temperature (Plate Temp Limit) automatically will be set to its maximum allowed value, to guarantee a short heat-up time. In some cases this is not desirable, but the set limitation of the hotplate temperature (Plate Limit, see section 8.1.2) should not be overwritten after changing the temperature of the medium (Probe Temp). This automatic function can be switched off here.

To switch the automatic setting of the hotplate temperature on or off, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Plate Auto Set".

**Adjustment range:** On, Off

**Default value:** On

### 12.2.8 Safety Auto Set

When this option is enabled, the device will automatically set the safety temperature 15°C above the programmed set temperature every time the setpoint is changed. If a Pt100 probe is connected, the set temperature is the set value of the heating medium (Probe Temp), otherwise it is the hotplate temperature (Plate Temp). Every time the Safety Temperature was set by the device automatically, a short message will be displayed to inform the user. After that, the user can readjust the Safety Temperature manually by using the setting knob.

To switch the automatic setting of the Safety temperature on or off, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Safety Auto Set".

**Adjustment range:** On, Off

**Default value:** On

### 12.2.9 Ask Volume

To optimize the heat-up time, whenever the device is switched on or a Pt100 probe is connected to the device, the user will be asked automatically for the volume of the heating medium. This automatic query can be switched off in this menu item. If this function is turned off, the volume of liquid can only be changed manually. For more information, see section 8.1.1-Setting up the Heating Volume.

To switch the automatic Ask Volume function on or off, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Ask Volume".

**Adjustment range:** On, Off

**Default value:** On

### 12.2.10 Motor Auto On

When this option is enabled, the motor will be automatically switched on every time the hotplate is switched on. However the motor will not be switched off if the hotplate is switched from On to Off. If this function is disabled, the motor must be switched on separately.

To switch the automatic Motor Auto On function on or off, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Motor Auto On".

**Adjustment range:** On, Off

**Default value:** Off

### 12.2.11 Safety Stir Time

After device malfunction (e.g. hotplate overheating) the hotplate switches automatically off. But the device continues stirring the liquid for a certain time, to prevent heat accumulation and helping cooling down the liquid. The user can define the time for continue stirring (Safety Stir Time) in this setting. For more information, see section 11.1-Safety Stir Function.

To adjust the Safety Stir Time, enter the Settings Menu (Section: 12.2) and choose the setting parameter "Safety Stir Time".

**Adjustment range:** 00:00:00 ... 01:00:00 h:m:s (off: 00:00:00 h:m:s)

**Default value:** 00:05:00

## 12.3 PID Menu

In the PID Menu you can change various setting parameters which affect the heating control. It is recommended that these setting parameters should only be changed by experienced users. This submenu can be reached over the Setup Menu as described in Section 12.1-Setup Menu and contains the following setting parameters:

- Thermal Resist. (Section: 12.3.1-Thermal Resistance)
- Container Type (Section: 12.3.2-Container Type)
- BACK

### 12.3.1 Thermal Resistance

This value reflects the energy drop from the liquid to ambient temperature.

This means that if there is a lot of energy needed (which comes from the hotplate) to hold the desired liquid temperature (Probe Temp), there is a low thermal resistance present. If for example water should be heated to 99°C at normal ambient conditions, this would need much more than only the double the energy to heat water to 50°C. If, however, a temperature close to the boiling temperature should be controlled, this value should be decreased to lower values to achieve a faster heat up time.

To adjust the Thermal Resistance of the heating medium enter the PID Menu (Section: 12.3) and choose the setting parameter "Thermal Resist."

**Adjustment range:** 50 ... 400

**Default value:** 380

#### **WARNING**

The value of the thermal resistance should only be readjusted from default when operating close to the boiling point of the liquid. Otherwise an extreme temperature overshoot can occur!

### 12.3.2 Container Type

This option allows to adapt/optimize the temperature regulation circuitry to the type of container (=vessel) being used.

To adjust the Container Type of the vessel, enter the PID Menu (Section: 12.3) and choose the setting parameter "Container Type".

**Adjustment range:** Glass, Aluminium, Stainless steel

**Default value:** Glass

## 12.4 Calibration Menu

The device allows a recalibrating of the Pt100 probe. This can be done by an easy to use 2-point calibration. The calibration submenu can be reached over the Setup Menu as described in Section 12.1-Setup Menu and contains the following setting parameters:

- Pt100A low Cal. (Section: 12.4.1-Calibration for Pt100A)
- Pt100A high Cal (Section: 12.4.1-Calibration for Pt100A)
- Reset Calibration (Section: 12.4.2-Reset Calibration)
- BACK

### 12.4.1 Calibration for Pt100A

The device allows you to recalibrate the Pt100 probe. To achieve a good calibration result, please operate carefully.

Please follow strictly to the prescribed calibration instructions. Prerequisite for a good result is a highly accurate thermometer:

- Enter the Calibration Menu as described in section 12.4.
- If already a calibration for the Pt100 probe exists, please reset it as described in chapter 12.4.2.
- Put the Pt100 probe together with the sensor of a high-precision temperature measuring device in ice water (both at least 5cm immersed) and wait at least 5 minutes until the temperature sensors has reached the temperature of the water.
- Select the menu item "Pt100A low Cal." from the Calibrate Menu.
- Now read the ice water temperature from the high-precision temperature measuring device. Set this temperature value in the bottom line of the display from the magnetic stirrer. In the upper line you will see the temperature value measured by the magnetic stirrer of the ice water (before calibration).
- The set value is entered automatically after a short time or can be confirmed manually by pressing the Timer key.
- Put the Pt100 probe together with the sensor of a high-precision temperature measuring device in boiling water (both at least 5cm immersed) and wait at least 5 minutes until the temperature sensors has reached the temperature of the water.
- Select the menu item "Pt100A high Cal." from the Calibrate Menu.
- Now read the temperature of the boiling water from the high-precision temperature measuring device. Set this temperature value in the bottom line of the display from the magnetic stirrer. In the upper line you will see the temperature value measured by the magnetic stirrer of the boiling water (before calibration).
- The set value is entered automatically after a short time or can be confirmed manually by pressing the Timer key.

To be sure that the new calibration is working properly, you have to control the result. The best way to do this is by heating a suitable medium at different temperatures (see section 8.1-Temperature Control with External Pt100 Probe) with the magnetic stirrer. Wait until the several temperatures are stable, before comparing the measured temperature of the magnetic stirrer with the temperature of a high-precision temperature measuring device.

### 12.4.2 Reset Calibration

To reset the user calibration and returning to the factory calibration of the Pt100 probe, enter the Calibration Menu (Section: 12.4) and choose the setting parameter "Reset Calibr."

**Adjustment range:** No, Pt100A (reset Pt100 probe)

**Default value:** No

## 12.5 Default Settings Menu

This submenu allows to restore the device to its default settings. All by the user adjusted parameters will be lost. The Default Settings submenu can be reached over the Select Menu as described in Section 12.1-Setup Menu and contains the following setting parameters:

- Default Settings (Section: 12.5.1-Default Settings)
- BACK

### 12.5.1 Default Settings

By selecting this function, the device can be reset to its factory settings. It will erase all user settings.

To reset the device to its default settings, enter the Default Settings Menu (Section: 12.5) and choose the setting parameter "Default Settings".

**Adjustment range:** No, Yes

**Default value:** No

## 13. Display Messages

### 13.1 Messages at Starting

At switching on the device displays different system information to inform the user about different settings:

Display	Description
Last Off Cond.: XXX	If the device is switched off automatically, the Off Condition "XXX" is displayed at the next switching on, to inform the user. See chapter 13.4-Off Conditions.
Temp. is Shown in Fahrenheit	The temperature unit is changed. The temperatures are shown in Fahrenheit instead of Celsius.
Device Uses a Pt100 User Cal.	The user has stored a user calibration for the Pt100 probe connector.
Ramp is Set to XXX°C/h	A temperature ramp was set for the temperature control.
Thermal Resist. Is set to: XXX	The Thermal Resistance has been changed in the PID menu.
Container Type: XXX	The Vessel Type has been changed in the PID menu.
Max. Plate Temp. Limited to: XXX°C	The max. adjustable hotplate temperature has been limited in the Settings menu to the value "XXX".
No Pt100 Probe Connected	At switching on, no Pt100 probe was found. The device is in the hotplate control mode.

## 13.2 Messages at Shutdown

At switching off the device displays different system information to inform the user about different settings:

Display	Description
Remaining Stir Time: XX:XX:XX	The device is in the Safety Stir Mode (see chapter 11.1-Safety Stir Function). This means the hotplate is disabled and the device stirs the duration of the adjusted Safety Stir Time before it turns off automatically.
Off Condition: XXX	If the device is switched off automatically, the Off Condition "XXX" is displayed to inform the user. See chapter 13.4-Off Conditions.
Hotplate is hot! Temp: XXX°C	The device was switched off, but the hotplate is still hotter than 60°C. You may get burned. This warning message is displayed until the hotplate is cooled below 60°C.

## 13.3 Information and Warnings

During Operation, the device displays different information to the user to simplify the handling of the device. Additionally, in the most cases a warning message is displayed to facilitate troubleshooting. A description of each message on the screen is described below:

Display	Description
Timer Expired	The Timer function was expired.
Multi Timer Expired	The Multitimer function was expired.
Multi Timer: No Steps Defined	The Multitimer function has been started, but there are no programmed steps.
Multi Timer: No Pt100 Probe	In at least one step in the Multitimer function, the Pt100 has been programmed but it is not connected to the device.
Plate and Motor Switched On	The hotplate and the motor were turned on automatically.
Probe DiffAlarm Has Triggered	The Differential Alarm was triggered. The temperature drop at the Pt100 probe was too quick. Maybe the Pt100 probe fell out of the vessel or the liquid has leaked.
OutOfLiq. Has Triggered	The OutOfLiquid Function was triggered. Maybe the external Pt100 probe is not immersed into the medium.
Probe Reached Safety Temp.	The measured temperature of the Pt100 probe was higher than the set Safety Temperature.
Set Probe Safety to: XXX°C	The safety temperature of the Pt100 probe was automatically set to the value "XXX"°C.
Probe Calibration Err.	The calibration of the Pt100 probe is damaged.
Probe Connected	The Pt100 probe was connected.
Probe Broken	The Pt100 probe is damaged or was disconnected during operation.
Contact Thermo. Connected	The Contact Thermometer was connected.
Contact Thermo. Disconnected	The Contact Thermometer is damaged or was disconnected during operation.
Plate Reached Over Temp.	The measured hotplate temperature was higher than the max. permissible temperature range.

Plate Reached Safety Temp.	The measured hotplate temperature was higher than the set Safety Temperature.
Set Plate Safety to: XXX°C	The safety temperature of the hotplate was automatically set to the value "XXX"°C.
Plate Calibration Err.	The calibration of the hotplate temperature sensor is damaged.
Plate Broken	The hotplate temperature sensor is damaged.
Plate Amplifier Shorted	The hotplate output stage is damaged.
Set Plate to Max: XXX°C	The hotplate temperature limit was automatically set to the max. possible value "XXX"°C.
Set Plate to Min: XXX°C	The hotplate temperature limit was automatically set to the min. possible value "XXX"°C.
Plate Limit Might Be To Low	Maybe the hotplate temperature limit is set too low to heat up the heating medium to the set temperature.
Front Internal Comm. Error	There was an error in the communication between the printed circuit boards inside the device.
Motor Internal Comm. Error	There was an error in the communication between the printed circuit boards inside the device.
Internal Over Temperature	The temperature inside the device has reached its maximum allowed temperature.
Internal Temp. Error	The internal temperature sensor is damaged.
Locked-Remote Control Mode	The control of the device via the control panel was disabled via RS485. Only the I/O key still works.

### 13.4 Off Conditions

Because of the many features and security mechanisms of the device, there are several possibilities why the device switches off itself automatically. To inform the user about the Off Condition and simplify locating and correcting problems, an according message appears in the display when the device switches off. The next time it is switched on, the last Off Condition is displayed again. If the device is controlled via the RS485 interface, the last Off Condition Code can be read with the RAC Command. A description of each Off Conditions on the screen and RS485 Off Condition Codes can be found below:

Display	RS485 Off Condition Code	Description
Switch Off	101	The device was switched off by the I/O key on the control panel.
Remote Off	102	The device was switched off by RS485 command.
Timer Expired	103	The Timer function was expired.
Multi Timer Exp.	104	The Multitimer function was expired.
Probe DiffAlarm	107	The Differential Alarm was triggered. The temperature drop at the Pt100 probe was too quick. Maybe the Pt100 probe fell out of the vessel or the liquid has leaked.

Probe OutOfLiq.	108	The OutOfLiquid Function was triggered. Maybe the external Pt100 probe is not immersed into the medium.
Probe Safety	109	The measured temperature of the Pt100 probe was higher than the set Safety Temperature.
Probe Broken	115	The Pt100 probe is damaged or was disconnected during operation.
ContactTherBroken	119	The Contact Thermometer is damaged or was disconnected during operation.
Plate OverTemp.	120	The measured hotplate temperature was higher than the max. permissible temperature range.
Plate Safety	122	The measured hotplate temperature was higher than the set Safety Temperature.
Plate Broken	127	The hotplate temperature sensor is damaged.
PlateAmpShort.	132	The hotplate output stage is damaged.
FrontIntComErr	136	There was an error in the communication between the printed circuit boards inside the device.
MotorIntComErr	137	There was an error in the communication between the printed circuit boards inside the device.
EEPROM Error	138	Faulty data was found in the EEPROM. The device tries to repair the data.
Intern Temp Err.	141	The temperature inside the device has reached its maximum allowed temperature or the internal temperature sensor is damaged.
MainsVoltageErr	142	The device was not turned off properly. Maybe the power cable was disconnected during operation or there was a mains power failure.
Watchdog Timer	144	The internal Watch Dog Timer of the microcontroller was triggered.

## 14. RS485 Interface

The RS485 interface allows remote access of all functions as well as readout and change of all system parameters (e.g. hotplate temperature, motor speed etc.).

Multiple devices may be connected parallel on a RS 232 port of a PC with an RS485 to RS232 adapter (PN: 61703-50) and connection cords for connection of different devices via RS485.

Therefore it is necessary to address the device. Every device has a so-called slave address between 1 and 255 (default address: 1). For parallel operation, each device needs its own and unique address, to avoid data clashes.

## 14.1 Settings of the RS485 Interface


The parameters of the interface (default settings: 9600, 8, N, 1):

<b>Baudrate:</b>	<ul style="list-style-type: none"> <li>• 1200 Baud</li> <li>• 2400 Baud</li> <li>• 4800 Baud</li> <li>• 9600 Baud (default setting)</li> </ul>
<b>Databits:</b>	<ul style="list-style-type: none"> <li>• 8 Bit</li> </ul>
<b>Parity:</b>	<ul style="list-style-type: none"> <li>• none</li> </ul>
<b>Stopbits:</b>	<ul style="list-style-type: none"> <li>• 1</li> </ul>
<b>Slave Address:</b>	<ul style="list-style-type: none"> <li>• 1 ... 255 (1 = default setting)</li> </ul>

## 14.2 Pin assignment of the RS485 Connector

	<b>NOTICE</b>	<p>Only connect original accessories to the RS485 jacks! Any other use may result in damage of the device or damage to third party devices!</p>
---	---------------	---

The serial port is a RS485 4-wire system with 2 differential inputs (A and B) and 2 differential outputs (Y and Z). The RS485 connector on the back of the device has the following pin assignment:

### 14.3 Connection between PC and Instrument


#### **NOTICE**

Connect the plug of the RS485 to RS232 adapter only to the RS485 jack of the device. Any other use may result in damage of the adapter or damage to third party devices!

For connecting the device with a PC a separate RS485 to RS232 adapter (PN: 61703-50) is necessary. You have to connect the RS485 plug of this adapter with one of the RS485 jacks on the rear of the device. The 9-pin Sub-D connector of the adapter is for the connection to the PC and has the following pin assignment.

Pin number	Description
2	transmit data TxD
3	receive data RxD
5	ground, GND

9-pin Sub-D connector (female) of the RS485 to RS232 adapter (view on socket)

For the extension of the cable you can use a standard 9-pin cable (pins not crossed, PN: 30275-51). When your PC has no RS232 port, you can use additionally an USB to RS232 converter (PN: 30244-01).

### 14.4 Addressing the Instrument

To use several devices on a bus-system, every device needs its own and unique address number. This address number is used to send the commands to the device. No address number must occur twice in the bus system to avoid data collision. The address numbers are possible between 1 and 255. The default setting is 1 for all devices. To change the address number, follow the instructions described in chapter 12.2.1-Slave Address.

## 14.5 RS485 Data Transfer

Every data transfer is started from the master (normally the PC). The master sends a command to the device (slave). After that the slave sends an answer back to the master, consisting of the repeated command and a handshake.

### Note:

The master must not send the next command until the complete answer from the slave is received (repeated command and the handshake)!

### 14.5.1 Format of an RS485 Command

Every command corresponds to the following format:

**ADR,CMDCODE,PARAMETERLIST<CR>**

Description:

- **ADR:** Slave Address of the desired device
- **CMDCODE:** Command-code
- **PARAMETERLIST:** 1 to 6 parameters separated by commas (see command table)
- **CR:** The command string must be terminated by Carriage/Return (ASCII code 13)

### 14.5.2 Format of an RS485 Handshake

When a slave device receives a command with its address it will answer to the master in this way:

- The complete command is send back
- The slave device sends a handshake corresponding to the following format:

**ADR,"HS",RETURNCODE,PARAMETERLIST<CR>**

Description:

- **ADR:** Slave Address of the Controller sending the handshake
- **RETURNCODE:** Error-code (see table below)
- **PARAMETERLIST:** 0 to 6 parameters (see command table), each parameter is separated by a comma ","
- **CR:** the handshake as any command, is terminated by ASCII-code 13 (CR)

Return Code	Explanation	Parameterlist
OK	command executed, no error	see command table
UC	unknown command	None
PA	wrong parameter number (too few or too many parameters specified)	None
NA	command is not allowed in actual operation mode	actual operation mode
PR	at least one parameter is out of range	None
PL	at least one parameter is too long	None
DF	unknown data format	None

## 14.6 RS485 Commands

Overview of all available RS485 Commands:

Command	Function	Parameter Description	Parameter Range	Example	Comment
<b>RTY</b>	Read Type and Version of device	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. name/type of device 2. Version number of software 3. On Off counts 4. Total Minutes Operation time	1  text number number number	1,RTY,1	
<b>PON</b>	Switch on Device	1. Security parameter 1234	1234	1,PON,1234	
<b>OFF</b>	Switch off Device	1. Security parameter 1234	1234	1,OFF,1234	
<b>WON</b>	Set Status (ON/OFF control) of - Motor - Plate	1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	0/1 0/1	1,WON,1,0	
<b>RON</b>	Read status of:  - Motor - Plate	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	1  0/1 0/1	1,RON,1	
<b>RAC</b>	Read actual	1. Dummy parameter to initiate transfer  -> Controller sends in handshake 1. actual motor speed in rpm 2. actual plate temperature in °C or °F 3. actual probe temperature in °C or °F (if connected, if not 'x' as return value) 4. actual safety Probe temperature in °C or °C (if connected, if not 'x' as return value) 5. last Off Condition Code	1  0...MaxSpeed* 0...MaxPlateTemp* 0...MaxProbeTemp*, x 0...MaxProbeTemp*, x 101...143** see: table Chapter 13.4.	1,RAC,1	Motor disabled in H30/30D. Device will send 'x' as parameter.  Safety probe temperature disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M36, MCS77, MCS78. Device will send 'x' as parameter.
<b>WSE</b>	Write actual set values of: - Motor - Hotplate (plate)  - external Pt100 sensor (probe)	1. Setpoint of motorspeed in rpm 2. Setpoint of hotplate temperature in °C or °F  3. Setpoint of probe temperature in °C or °F	0, 60...MaxSpeed* <b>Without Pt100:</b> 0...MaxPlateTemp* <b>With Pt100:</b> SetProbeTemp+10°C... MaxPlateTemp* 0...MaxProbeTemp*	1,WSE,800,34 0,60	Motor disabled in H30/30D. Device will ignore the value.
<b>RSE</b>	Read actual set values of:  - Motor - Hotplate (plate)  - external Pt100 sensor (probe)	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Setpoint of motorspeed in Rpm 2. Setpoint of hotplate temperature in °C or °F 3. Setpoint of probe temperature in °C or °F	1  0...MaxSpeed* 0...MaxPlateTemp* 0...MaxProbeTemp*	1,RSE, 1	Motor disabled in H30/30D. Device will send 'x' as parameter.
<b>WTR</b>	Set / Write:  - Timer  - Ramp  - Safety temperature	1. Timer value in seconds. Set this value to 0 to disable the timer. 2. setting of the ramp in °C/h or °F/h (a value of 450°C/h disables the ramp) 3. Safety temperature in °C	0...MaxTimer*  1...450  <b>Without Pt100:</b> -PlateTemp+1... MaxPlateTemp*+25°C <b>With Pt100:</b> -ProbeTemp + 1... MaxProbeTemp*+25°C	1,WTR,600,45 0,160	Ramp disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will ignore the value.

<b>RTR</b>	Read settings of:  - Timer  - Ramp  - Safety temperature	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. actual setting of the timer in seconds. If the timer value is 0, the timer is disabled. 2. setting of the ramp in °C/h (a value of 450 signals that the ramp is disabled) 3. Safety temperature in °C	1  0...MaxTimer*  1...450, x  <b>Without Pt100:</b> -0...MaxPlateTemp*+ 25°C <b>With Pt100:</b> -0...MaxProbeTemp*+ 25°C	1,RTR,1	Ramp disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will send 'x' as parameter.
<b>WMS</b>	Set / Write Multitimer Step Parameter: -Step number -Timer  -Hotplate temp.  -Probe temp.  -Ramp  -Motor speed	1. Step number 2. Timer value in seconds. Set this value to 0 to disable this timer step. Set it to -1, -2 or -3 for plate reached, probe reached or motor reached setpoint 3. Setpoint of hotplate temperature in °C or °F  4. Setpoint of probe temperature in °C or °F 5. Setting of the ramp in °C/h or °F/h (a value of 450°C/h disables the ramp) 6. Setpoint of motorspeed in rpm	1...5 -3...MaxTimer*  <b>Without Pt100:</b> 0...MaxPlateTemp* <b>With Pt100:</b> SetProbeTemp+10°C... MaxPlateTemp* 0...MaxProbeTemp*  1...450  0... MaxSpeed *	1,WSM,2,600, 200,100,10,2 00	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>RMS</b>	Read Multitimer Step Parameter:  -Step number -Timer  -Hotplate temp.  -Probe temp.  -Ramp  -Motor speed	1. Step number  -> Controller sends in handshake: 1. Step number 2. Timer value in seconds. If this value is set to 0 this timer step is disabled. If it is set to -1, -2 or -3, then plate reached, probe reached or motor reached setpoint is activated. 3. Setpoint of hotplate temperature in °C or °F  4. Setpoint of probe temperature in °C or °F 5. Setting of the ramp in °C/h (a value of 450 signals that the ramp is disabled) 6. Setpoint of motor speed in rpm	1...5  1...5 -3...MaxTimer*  <b>Without Pt100:</b> 0...MaxPlateTemp* <b>With Pt100:</b> SetProbeTemp+10°C... MaxPlateTemp* 0...MaxProbeTemp*  1...450  0... MaxSpeed *	1,RMS,2	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>WMO</b>	Set / Write Multitimer Option Parameter: -Number of cycles  -Expired behavior	1. Number of Multitimer cycles. Set this value to 0 to repeat the steps endless. 2. Expired behavior: -0=Hold settings -1=Plate off -2=Plate and motor off -3=Switch device off	0...999  0...3	1,WMO,10,0	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>RMO</b>	Read Multitimer Option Parameter  -Number of cycles	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Number of Multitimer cycles. If this value is set to 0 the steps are	1  0...999	1,RMO,1	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36

	-Expired behavior	repeated endlessly. 2. Expired behavior: -0=Hold settings -1=Plate off -2=Plate and motor off -3=Switch device off	0...3		
<b>WT2</b>	Set Status (ON/OFF control) of Multitimer	1. On/Off (0-> Off, 1-> On)	0/1	1,WT2,1	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>RT2</b>	Read status of Multitimer:  -Status  -Actual cycles  -Actual step -Actual step time  -Running time	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Multitimer On/Off (0-> Off, 1-> On) 2. Actual number of Multitimer cycles 3. Actual Multitimer Step 4. Remaining Time of the actual Step 5. Running time	1  0/1  0...999  1...5 0...MaxTimer*  0...MaxTimer*	1,RT2,1	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>WVO</b>	Write/Set volume	1. Volume in ml	100...9900	1,WVO,1000	
<b>RVO</b>	Read volume	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. programmed volume in ml	1  100...9900	1,RVO,1	
<b>WSM</b>	Activate serial mode (→ enable/disable setpoint control via frontpanel)	1. Set this value to 0 to allow control via the frontpanel or set this value to 1 to disable control via the frontpanel. (only used for devices with motor)	0 / 1		
<b>WTU</b>	Set / Write Temperature Units	1. Type of Units (0->Celsius / 1->Fahrenheit)	0 / 1	1,WTU,0	
<b>RTU</b>	Read Temperature Units	1. Dummy parameter to initiate transfer -> Controller sends in handshake: 1. Type of Units	1  0=Celsius / 1=Fahrenheit	1,RTU,1	
<b>WSU</b>	Set / Write Safety Auto Set function	1. Set this value to 0 to disable Safety Auto Set or set this value to 1 to enable Safety Auto Set. (0->Off / 1->On)	0 / 1	1,WSU,0	Safety Auto Set disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will ignore the value.
<b>RSU</b>	Read Safety Auto Set function	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Safety Auto Set On/Off (0-> Off, 1-> On)	1  0 / 1		Safety Auto Set disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will send 'x' as parameter.
<b>WSD</b>	Write setup data: - Plate Temp Limit - Safety Stir Time - Ask for Volume at power up - Differential Alarm Sensitivity - Out of Liquid Sensitivity - Thermal Resistance	1. Plate Temp Limit in °C 2. Safety Stir Time in seconds 3. Ask for Volume at power up (0-> Off, 1-> On) 4. Differential Alarm Sensitivity in % 5. Out of Liquid Sensitivity in % 6. Thermal Resistance	50...MaxPlateTemp* 0...3600 0 / 1  1...100  0...100 (0=Off)  50...400	1,WSD,360,600,1,90,20,380	Command completely disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will ignore the value.
<b>RSD</b>	Read setup data:  - Plate Temp Limit - Safety Stir Time - Ask for Volume at power up	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Plate Temp Limit in °C 2. Safety Stir Time in seconds 3. Ask for Volume at power up (0-> Off, 1-> On)	1  50...MaxPlateTemp* 0...3600 0 / 1	1,RSD,1	Command completely disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will send 'x' as parameters.

	- Differential Alarm Sensitivity - Out of Liquid Sensitivity - Thermal Resistance	4. Differential Alarm Sensitivity in % 5. Out of Liquid Sensitivity in % 6. Thermal Resistance	1...100  0...100 (0=Off)  50...400		
<b>RCO</b>	Read connector status	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Connector probe status 0 - Nothing connected 1 - PT100 probe connected 2 - PT100 dummy connected 3 - Contactthermometer mode 2. Connector Safety probe status 0 - Nothing connected 1 - PT100 probe connected 2 - PT100 dummy connected 3 - Contactthermometer mode	1  0...3  0...3	1,RCO,1	Connector Safety probe is disabled in M21, M22, M23 KM16.4, KM16.7, M36, H30/30D, MCS77 and MCS78. Device will send 'x' as parameters.
<b>RSS</b>	Read system status	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Device On/Standby/Safety Stir  2. Remaining Safety Stir Time in seconds	0=Standby / 1=On / 2=Safety Stir Mode 0...3600	1,RSS,1	
<b>WSA</b>	Set RS485 slave address	1. New slave address of controller	0...255	1,WSA,3	
<b>WBD</b>	Set RS485 Baudrate	1. Baudrate	0=1200baud / 1=2400baud / 2=4800baud / 3=9600baud	1,WBD,2	
<b>RST</b>	Reset the device to default settings	1. Security parameter 1234	1234	1,RST,1234	

**Note:**

In the Fahrenheit mode all values are shown in Fahrenheit and not in Celsius!

\*The RS485 commands are applicable to different device types with different parameter ranges. In the following table you can find all device specific parameters of all device types:

Device Type	MaxPlateTemp	MaxProbeTemp	MaxSpeed	MaxTimer
M 21	350°C	250°C	1600rpm	59940
M 22	380°C	250°C	1600rpm	59940
M 23	500°C	250°C	1600rpm	59940
KM 16.4	450°C	300°C	1100rpm	59940
KM 16.7	450°C	300°C	1100rpm	59940
H 30/30D	380°C	250°C	disabled	59940
M 26G2	360°C	250°C	1600rpm	86400
M 36	500°C	250°C	1100rpm	86400
MCS 77	330°C	250°C	1600rpm	86400
MCS 78	440°C	250°C	1600rpm	86400

\*\*The appropriate assignment of the "RS485 Off Condition Codes" to the corresponding Off Conditions can be found in the table in chapter 13.4-Off Conditions.

## 14.7 RS485 Examples

### 14.7.1 Example 1 – Setting up a Value

Programming of a MCS 77 device with the following parameters:

- Hotplate temperature: 300°C (corresponds to the plate limit)
- Probe temperature: 50°C
- Motor speed: 500rpm

The device has the standard slave address 1 and a Pt100 probe is connected:

Following commands must be send to the device. The answers of the MCS 77 are marked *italic*, they must not be send to the device:

<b>1,PON,1234&lt;CR&gt;</b>	;Switch on the device.
<b>1,PON,1234&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1,HS,OK&lt;CR&gt;</b>	;Handshake sent from the device, command accepted.
<b>1,WSM,1&lt;CR&gt;</b>	;Deactivating the front panel only controlling via RS485 (optional).
<b>1,WSM,1&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1,HS,OK&lt;CR&gt;</b>	;Handshake sent from the device, command accepted.
<b>1,WSE,500,300,50&lt;CR&gt;</b>	;Set values (motor speed, hotplate and probe temperature).
<b>1,WSE,500,300,50&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1, HS,OK&lt;CR&gt;</b>	;Handshake sent from the device, command accepted.
<b>1,WON,1,1&lt;CR&gt;</b>	;Switch on motor and hotplate.
<b>1,WON,1,1&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1, HS,OK&lt;CR&gt;</b>	;Handshake sent from the device, command accepted.

## 14.7.2 Example 2 – Reading a Value

Reading of the actual values of a MCS 77 device:

- Temperature units
- Hotplate temperature
- Probe temperature
- Motor speed

The device has the standard Slave Address 1 and a Pt100 probe is connected:

Following commands must be sent to the device. The answers of the MCS 77 are marked *italic*, they must not be sent to the device:

<b>1,PON,1234&lt;CR&gt;</b>	;Switch on the device.
<b>1,PON,1234&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1,HS,OK&lt;CR&gt;</b>	;Handshake sent from the device, command accepted.
<b>1,RTU,1&lt;CR&gt;</b>	;Reading temperature unit.
<b>1,RTU,1&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1,HS,OK,0&lt;CR&gt;</b>	;Handshake sent from the unit, command accepted. (temperature unit = Celsius)
<b>1,RAC,1&lt;CR&gt;</b>	;Reading actual values.
<b>1,RAC,1&lt;CR&gt;</b>	;Command is sent back repeatedly from the device.
<b>1,HS,OK,480,180,50,x,101&lt;CR&gt;</b>	;Handshake sent from the unit, command accepted. (motor speed = 480rpm, hotplate temperature = 180°C, probe temperature = 50°C)

The actual settings can be read out of the different handshakes from the slave device:

- The temperatures are shown in Celsius (command RTU, parameter 1 = 0)
- The actual motor speed is 480rpm (command RAC, parameter 1 = 480)
- The actual hotplate temperature is 180°C (command RAC, parameter 2 = 180)
- The actual probe temperature is 50°C (command RAC, Parameter 3 = 50)

## 15. Maintenance and Cleaning

### **WARNING**

In case of malfunction do not attempt to repair the device.  
There are no user-serviceable parts in this instrument.

The outer casting is mostly coated aluminum, the hotplate consist of anodized aluminum (MCS77) or Ceran (MCS78). Together with a chemically resistant splash-proof membrane key-pad and therefore easily cleaned with warm water and any proprietary liquid laboratory detergent. Do not use steel wool or any similar plastic wool sponge or any aggressive cleaning agents to clean the device.

Crust, limestone and water stains may be cleaned with the usual household detergents. Please remove remaining detergent thoroughly with a damp cloth or sponge, since some of the detergents react caustic with higher temperatures. If soils are already encrusted use a razorblade scratcher. Plastics, kitchen foil as well as stuff containing sugar must be removed immediately.

Mechanical supports are: Razorblade scratchers, non-scratching sponges and non-scratching steel wool.

The device should only be opened and repaired by authorized service personnel. Any work on the electronics in the device should only be carried out by knowledgeable, trained personnel. Any attempt by the user to repair the device will immediately render the guarantee null and void. Please contact your local distributor in the event of a problem.

## 16. Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's election, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement device will be warranted for the remainder of the original twenty-four (24) month period or ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein.

The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damage of any kind resulting from use or possession of the product.

### **WARNING**

The onus is on the user to find out whether the device is suitable for his application. If in doubt clarify this with your dealer or with the manufacturer directly.

## 17. Dismantling, Transport and Storage

### 17.1 Dismantling

- Switch off the unit by pressing the I/O key.
- Disconnect the unit from the mains.
- Now you may remove the instrument from the working area.

### 17.2 Transport and Storage

#### Prior to transport:

- Switch the instrument off and unplug the power supply.
- Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.
- Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).
- Do not subject the instrument to mechanical shocks or vibration during transporting it.
- In case you do not use the original packaging please mark the box with the following notes:
  - Glass symbol (handle with care, fragile)
  - Umbrella (keep dry)
  - Content (list of content)

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, please sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

WEEE No. DE 54398924

## 18. Technical Data

Technical Data MCS 77 / 78	
Operating the device:	<ul style="list-style-type: none"><li>• Control panel</li><li>• RS485</li></ul>
Display:	<ul style="list-style-type: none"><li>• backlit graphical LCD display</li></ul>
Hotplate:	<p><b>MCS 77:</b></p> <ul style="list-style-type: none"><li>• Anodized aluminum</li><li>• Round, diameter 140mm</li><li>• Power: 500 Watt</li><li>• Adjustable temperature range: RT ... 330°C</li></ul> <p><b>MCS 78:</b></p> <ul style="list-style-type: none"><li>• Ceran</li><li>• Round, diameter 130mm</li><li>• Power: 600 Watt</li><li>• Adjustable temperature range: RT ... 440°C</li></ul> <ul style="list-style-type: none"><li>• Adjustable hotplate temperature (without Pt100)</li><li>• Adjustable max. hotplate temperature (with Pt100)</li></ul>

	<ul style="list-style-type: none"> <li>• Microprocessor controlled temperature control</li> <li>• Setting accuracy: 1°C</li> </ul>
<b>Pt100 Connector:</b>	<ul style="list-style-type: none"> <li>• Connector for an external Pt100 probe</li> <li>• Adjustable Pt100 probe temperature</li> <li>• Microprocessor controlled temperature control</li> <li>• Adjustable temperature range: RT ... 250°C</li> <li>• Setting accuracy: 0,2°C</li> <li>• Programmable temperature ramp rate: 1°C/hour to 450°C/hour</li> <li>• Possibility of Pt100 probe calibration</li> </ul>
<b>Motor:</b>	<ul style="list-style-type: none"> <li>• Shaded-pole motor</li> <li>• Power: 25 Watt</li> <li>• Adjustable speed</li> <li>• Microprocessor controlled speed control</li> <li>• Adjustable speed range: 60 ... 1600min<sup>-1</sup></li> <li>• Setting accuracy: 10min<sup>-1</sup></li> <li>• Speed control stability: +/- 20min<sup>-1</sup></li> <li>• Soft motor start</li> <li>• Adjustable safety stir function</li> <li>• Recommended stirring bar: -Length: 60mm (min: 30 ... max: 70mm), -Diameter: 8mm (min: 5 ... max: 12mm)</li> </ul>
<b>Timer:</b>	<ul style="list-style-type: none"> <li>• Adjustable switch off time</li> <li>• Adjustable time range: 24:00:00 h:m:s</li> <li>• Setting accuracy: 1 min</li> <li>• Multitimer function</li> </ul>
<b>Safety functions:</b>	<ul style="list-style-type: none"> <li>• Programmable hotplate safety temperature (without Pt100)</li> <li>• Programmable probe safety function (with Pt100)</li> <li>• Detects and protects against the following hazardous situations: <ul style="list-style-type: none"> <li>○ Hotplate failure (hotplate exceeds max. allowed temperature)</li> <li>○ Pt100 probe failure (disconnection or break of Pt100 probe)</li> <li>○ Heating stage failure</li> <li>○ Extreme rise or fall of temperature of the external Pt100 probe, e.g. drop out of the Pt100 probe or breakage of the glass (Differential Alarm)</li> <li>○ Liquid detection of the Pt100 probe at switch on of the heating function (Out of Liquid)</li> <li>○ Internal overheating of the device</li> </ul> </li> </ul>
<b>Electrical power requirements:</b>	<ul style="list-style-type: none"> <li>• 115/230 Volts, 50/60 Hz (see type label on rear panel)</li> </ul>
<b>Power consumption:</b>	<p><b>MCS 77:</b></p> <ul style="list-style-type: none"> <li>• 545 Watt</li> </ul> <p><b>MCS 78:</b></p> <ul style="list-style-type: none"> <li>• 645 Watt</li> </ul>
<b>Permissible ambient temperature:</b>	<ul style="list-style-type: none"> <li>• 5°C ... 40°C</li> </ul>
<b>Permissible humidity:</b>	<ul style="list-style-type: none"> <li>• 80% RH</li> </ul>
<b>Safety class:</b>	<ul style="list-style-type: none"> <li>• IP32 (according DIN 40050)</li> </ul>
<b>Dimensions:</b>	<ul style="list-style-type: none"> <li>• 180 x 115 x 250 mm (B x H x T)</li> </ul>
<b>Weight:</b>	<ul style="list-style-type: none"> <li>• 3,9 kg</li> </ul>

## 19. Repair Return Form

<b>IMPORTANT:</b>
-Calibration Service is automatically performed for instruments that require adjustments
-We contact you by sending a repair estimate.
-Repairs will only be carried out after your release.

<b>CONTACT / USER INFORMATION:</b>	
Contact:	Phone No.:
Fax No.:	E-Mail:
Billing:	Shipping:
Company:	Company:
Address:	Address

<b>INSTRUMENT INFORMATION:</b>	
Model:	Serial-No.:
Please describe all problems / malfunctions:	

<b>OPERATION CONDITIONS</b> (please fill in if applicable):		
Ambient temp. in °C:	Humidity:	Speed:
Load:	Volume:	Viscosity:
Temperature in °C	Sample temperature in °C:	Operating time:
Sample Description*:		

**\*NOTE:** If the instrument was exposed to hazardous material, it must be decontaminated **before** returning it to Ingenieurbüro CAT, M. Zipperer GmbH and the MSDS for hazardous material must be included with the instrument.

<b>RETURN SHIPPING:</b>			
UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>

\*\*Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means:

<b>PACKAGING INSTRUCTIONS TO RETURN AN INSTRUMENT FOR REPAIR:</b>
-Remove all accessories (e.g. homogenizer tools, stirring paddles, stands, clamps, cables) from the instrument
-Clean excess testing material off the instrument/accessory
-Include MSDS sheets for all hazardous materials used with this instrument
-Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
-DO NOT send accessories unless there is a problem with them. When sending back, wrap each item to avoid contact with the instrument.
-Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions. Example: "Handle with care" or "Fragile- Delicate Instrument" and send to:
<b>Ingenieurbüro CAT</b> <b>M. Zipperer GmbH</b> <b>Service Department</b> <b>Etzenbach 16</b> <b>D-79219 Staufen</b>

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Command	Function	Parameter Description	Parameter Range	Example	Comment
<b>RTY</b>	Read Type and Version of device	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. name/type of device 2. Version number of software 3. On Off counts 4. Total Minutes Operation time	1  text number number number	1,RTY,1	
<b>PON</b>	Switch on Device	1. Security parameter 1234	1234	1,PON,1234	
<b>OFF</b>	Switch off Device	1. Security parameter 1234	1234	1,OFF,1234	
<b>WON</b>	Set Status (ON/OFF control) of - Motor - Plate	1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	0/1 0/1	1,WON,1,0	
<b>RON</b>	Read status of:  - Motor - Plate	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Motor On/Off (0-> Off, 1-> On) 2. Plate On/Off (0-> Off, 1-> On)	1  0/1 0/1	1,RON,1	
<b>RAC</b>	Read actual	1. Dummy parameter to initiate transfer  -> Controller sends in handshake 1. actual motor speed in rpm 2. actual plate temperature in °C or °F 3. actual probe temperature in °C or °F (if connected, if not 'x' as return value) 4. actual safety Probe temperature in °C or °C (if connected, if not 'x' as return value) 5. last Off Condition Code	1  0...MaxSpeed* 0...MaxPlateTemp*  0...MaxProbeTemp*, x  0...MaxProbeTemp*, x  101...143** see: table Chapter 31.4.	1,RAC,1	Motor disabled in H30/30D. Device will send 'x' as parameter.  Safety probe temperature disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M36, MCS77, MCS78. Device will send 'x' as parameter.
<b>WSE</b>	Write actual set values of: - Motor - Hotplate (plate)  - external Pt100 sensor (probe)	1. Setpoint of motorspeed in rpm 2. Setpoint of hotplate temperature in °C or °F  3. Setpoint of probe temperature in °C or °F	0, 60...MaxSpeed* <b>Without Pt100:</b> 0...MaxPlateTemp* <b>With Pt100:</b> SetProbeTemp+10°C... MaxPlateTemp* 0...MaxProbeTemp*	1,WSE,800,34 0,60	Motor disabled in H30/30D. Device will ignore the value.
<b>RSE</b>	Read actual set values of:  - Motor - Hotplate (plate)  - external Pt100 sensor (probe)	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Setpoint of motorspeed in Rpm 2. Setpoint of hotplate temperature in °C or °F 3. Setpoint of probe temperature in °C or °F	1  0...MaxSpeed* 0...MaxPlateTemp*  0...MaxProbeTemp*	1,RSE, 1	Motor disabled in H30/30D. Device will send 'x' as parameter.
<b>WTR</b>	Set / Write:  - Timer  - Ramp  - Safety temperature	1. Timer value in seconds. Set this value to 0 to disable the timer. 2. setting of the ramp in °C/h or °F/h (a value of 450°C/h disables the ramp) 3. Safety temperature in °C	0...MaxTimer*  1...450  <b>Without Pt100:</b> -PlateTemp+1... MaxPlateTemp*+25°C <b>With Pt100:</b> -ProbeTemp + 1... MaxProbeTemp*+25°C	1,WTR,600,45 0,160	Ramp disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will ignore the value.

<b>RTR</b>	Read settings of:  - Timer  - Ramp  - Safety temperature	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. actual setting of the timer in seconds. If the timer value is 0, the timer is disabled. 2. setting of the ramp in °C/h (a value of 450 signals that the ramp is disabled) 3. Safety temperature in °C	1  0...MaxTimer*  1...450, x  <b>Without Pt100:</b> -0...MaxPlateTemp*+ 25°C <b>With Pt100:</b> -0...MaxProbeTemp*+ 25°C	1,RTR,1	Ramp disabled in M21, M22, M23 KM16.4, KM16.7 and H30/30D. Device will send 'x' as parameter.
<b>WMS</b>	Set / Write Multitimer Step Parameter: -Step number -Timer  -Hotplate temp.  -Probe temp.  -Ramp  -Motor speed	1. Step number 2. Timer value in seconds. Set this value to 0 to disable this timer step. Set it to -1, -2 or -3 for plate reached, probe reached or motor reached setpoint 3. Setpoint of hotplate temperature in °C or °F  4. Setpoint of probe temperature in °C or °F 5. Setting of the ramp in °C/h or °F/h (a value of 450°C/h disables the ramp) 6. Setpoint of motorspeed in rpm	1...5 -3...MaxTimer*  <b>Without Pt100:</b> 0...MaxPlateTemp* <b>With Pt100:</b> SetProbeTemp+10°C... MaxPlateTemp* 0...MaxProbeTemp*  1...450  0... MaxSpeed *	1,WSM,2,600,200,100,10,200	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>RMS</b>	Read Multitimer Step Parameter:  -Step number -Timer  -Hotplate temp.  -Probe temp.  -Ramp  -Motor speed	1. Step number  -> Controller sends in handshake: 1. Step number 2. Timer value in seconds. If this value is set to 0 this timer step is disabled. If it is set to -1, -2 or -3, then plate reached, probe reached or motor reached setpoint is activated. 3. Setpoint of hotplate temperature in °C or °F  4. Setpoint of probe temperature in °C or °F 5. Setting of the ramp in °C/h (a value of 450 signals that the ramp is disabled) 6. Setpoint of motor speed in rpm	1...5  1...5 -3...MaxTimer*  <b>Without Pt100:</b> 0...MaxPlateTemp* <b>With Pt100:</b> SetProbeTemp+10°C... MaxPlateTemp* 0...MaxProbeTemp*  1...450  0... MaxSpeed *	1,RMS,2	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>WMO</b>	Set / Write Multitimer Option Parameter: -Number of cycles  -Expired behavior	1. Number of Multitimer cycles. Set this value to 0 to repeat the steps endless. 2. Expired behavior: -0=Hold settings -1=Plate off -2=Plate and motor off -3=Switch device off	0...999  0...3	1,WMO,10,0	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36
<b>RMO</b>	Read Multitimer Option Parameter  -Number of cycles	1. Dummy parameter to initiate transfer  -> Controller sends in handshake: 1. Number of Multitimer cycles. If this value is set to 0 the steps are	1  0...999	1,RMO,1	Command disabled in M21, M22, M23, KM16.4, KM16.7, H30/30D, M26G2, M36

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**EC - Declaration of conformity in accordance with the EEC machine directive  
2006/42/EG appendix II 1A**

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We declare under our sole responsibility that this product (see table) corresponds to the regulations and conforms with the standards or standardized documents:

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The agreement with further valid guidelines/regulations following for the product is explained:

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Applied one harmonized standards in particular:

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In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

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